

## ADV003 Adenovirus Primer Set Instructions for Use

<b>Product Overview</b>	
Name	Adenovirus Primer Set
Part Number	ADV003
Intended Use	The Adenovirus Primer Set is an analyte specific reagent (ASR) for the detection of Adenovirus species A, B, C, D, E, F and G DNA by polymerase chain reaction (PCR).
Regulatory Status	Analyte Specific Reagent. Analytical and performance characteristics are not established.
Number of vials	1
Fill Volume	600µL
Storage Temperature	-20°C or below
Number of acceptable of freeze thaws	24

<b>Product Information</b>	
Analyte Target	Adenovirus
Fluorophore	FAM
Sequence Target	Hexon
Organisms determined to not have any binding affinity via <i>in silico</i> analysis	BK Virus, CMV, EBV, JC Virus, HHV-6, VZV, and HSV
Handling	Always handle samples and product components with caution to avoid contamination. Protect product from prolonged exposure to light. Thaw completely, mix by pipetting or gentle vortexing and centrifuge briefly before use.
Precautions	Discard sample and product waste according to your local safety regulations.
Limitations	The product is for professional use only. Each laboratory must establish its own standard operating procedures and instructions.