

## BAT001 *Bacillus atrophaeus* Primer Set Instructions for Use

<b>Product Overview</b>	
Name	Bacillus atrophaeus Primer Set
Part Number	BAT001
Intended Use	The Bacillus atrophaeus Primer Set is an analyte specific reagent (ASR) for the detection of <i>Bacillus atrophaeus</i> DNA by polymerase chain reaction (PCR).
Regulatory Status	Analyte Specific Reagent. Analytical and performance characteristics are not established.
Number of vials	1
Fill Volume	192µL
Storage Temperature	-20°C or below
Number of acceptable of freeze thaws	24

<b>Product Information</b>	
Analyte Target	<i>Bacillus atrophaeus</i>
Fluorophore	HEX
Sequence Target	<i>atpD</i>
Handling	Always handle samples and product components with caution to avoid contamination. Protect product from prolonged exposure to light. Thaw completely, mix by pipetting or gentle vortexing and centrifuge briefly before use.
Precautions	Discard sample and product waste according to your local safety regulations.
Limitations	The product is for professional use only. Each laboratory must establish its own standard operating procedures and instructions.