

DIL100 Specimen Diluent Instructions for Use

Product Overview	
Name	Specimen Diluent
Part Number	DIL100
Intended Use	Specimen Diluent is a General Purpose Reagent (GPR) and is intended for sample preparation. The solution is free of substances that can inhibit PCR and can be used in combination with many sample types.
Regulatory Status	General Purpose Reagent. For Laboratory Use.
Number of vials	100
Fill Volume	975µL
Storage Temperature	15°C to 25°C

Product Information	
Component Information	10 mM Tris, 0.1 mM EDTA pH 8.0
Handling	Always handle samples and product components with caution to avoid contamination. Vials are provided without any labeling. Laboratory specific labeling should be added prior to use.
Precautions	Discard sample and product waste according to your local safety regulations.
Limitations	The product is for professional use only. Each laboratory must establish its own standard operating procedures and instructions.