

## PEV002 Parechovirus Primer Set Instructions for Use

<b>Product Overview</b>	
Name	Parechovirus Primer Set
Part Number	PEV002
Intended Use	The Parechovirus Primer Set is an analyte specific reagent (ASR) for the detection of Parechovirus RNA by polymerase chain reaction (PCR).
Regulatory Status	Analyte Specific Reagent. Analytical and performance characteristics are not established.
Number of vials	1
Fill Volume	600µL
Storage Temperature	-20°C or below
Number of acceptable of freeze thaws	24

<b>Product Information</b>	
Analyte Target	Parechovirus
Fluorophore	FAM
Sequence Target	UTR
Organisms determined to not have any binding affinity	Enterovirus, Influenza A H1N1, Influenza A H1N1 pdm09, Influenza A H3N2, Influenza B, SARS-CoV-2, Rhinovirus 1A, RSV
Handling	Always handle samples and product components with caution to avoid contamination. Protect product from prolonged exposure to light. Thaw completely, mix by pipetting or gentle vortexing and centrifuge briefly before use.
Precautions	Discard sample and product waste according to your local safety regulations.
Limitations	The product is for professional use only. Each laboratory must establish its own standard operating procedures and instructions.